



*Zebra*



*Zebra*

# Pediatric Exam 1 Adults



## PT20 - S Archimed Zebra Pediatrik Muayene Masası

Archimed Pediatrik Muayene Masası Zebra figürlü olup, çocuklar için renkli ve ilgi uyandıran polikliniklere tasarlanmıştır. Üst minder döşemesi, sırt, oturak ve ayak uzatma olarak 3 kısımdan oluşmaktadır. Sırt kısmı mekanik kademe mekanizma ile kolayca ayarlanabilir. 2 adet çekmece ile konforlu bir kullanım sağlar. Tekerlekli olduğundan taşıma kolaylığı sunar. Tek bir pedal ile masayı sabit ayak üzerine almak mümkündür. Üründe kullanılan boya ve döşeme antibakteriyel özelliktedir.



## PT20 -S Archimed Zebra Pediatric Examination Table

Archimed Pediatric Examination Table is designed as Zebra figure which is attract attention and colored to design for polyclinics.  
Upholstery consist of 3 main parts backrest, seat and extension part under the upholstery.  
Backrest can be adjusted easily with mechanically stepped mechanism.  
The table provides comfortable usage with 2 pieces drawer.  
The table provides easy move by wheels.  
The table can be fixed on the ground with a single pedal.  
The paint which use on table and upholstery has antibacterial feature.



Sirt hareketi  
Back rest movement

Uzama Parçası  
Extension Peace



Çekmece / Drawer



## Teknik Özellikler Technical Specifications

Komple Uzunluk - Full Length	1780 mm
En - Width	620 mm
Yükseklik Oturak - Seat Height	900 mm
Oturak Uzunluk - Seat Length	1140 mm
Ayak Uzatma - Foot Extension	370 mm
Sünger Dansite - Foam Density	28
Sırt Açısı - Back Angle	30°
Ayak Uzatma En - Foot Extension Width	400 mm
Ürün Ağırlığı - Product Weight	89 kg
Teker Çapı - Wheel Diameter	75 mm
Maksimum Yükseklik - Max. Height	150 mm



Arşimed Medikal Kozm. San. Tic Ltd. Şti.

Mimar Sinan Mah. Bosna Cd. No:5 / 3 Çekmeköy - İstanbul

☎ : +90 216 567 8 567 (pbx) 📠 : +90 216 567 8 367 (pbx)

🌐 : www.arsimed.com.tr 📧 : info@arsimed.com.tr



CE FDA EN 60601-1 EN 60601-2



diabetic exam tables