







PT30 - S Archimed Deve Pediatrik Muayene Masası



Archimed Pediatrik Muayene Masası Deve figürlü olup, çocuklar için renkli ve ilgi uyandıran polikliniklere tasarlanmıştır. Üst minder döşemesi, sırt, oturak ve ayak uzatma olarak 3 kısımdan oluşmaktadır. Sırt kısmı mekanik kademe mekanizma ile kolayca ayarlanabilir. 2 adet çekmece ile konforlu bir kullanım sağlar. Tekerlekli olduğundan taşıma kolaylığı sunar. Tek bir pedal ile masayı sabit ayak üzerine almak mümkündür. Üründe kullanılan boya ve döşeme antibakteriyel özelliindedir.



PT30 - S Archimed Camel Pediatric Examination Table

Archimed Pediatric Examination Table is designed as Camel figure which is attract attention and colored to design for polyclinics.

Upholstery consist of 3 main parts backrest, seat and extension part under the upholstery.

Backrest can be adjusted easily with mechanically stepped mechanism.

The table provides comfortable usage with 2 pieces drawer.

The table provides easy move by wheels.

The table can be fixed on the ground with a single pedal.

The paint which use on table and upholstery has antibacterial feature.



Çekmece / Drawer



Sirt hareketi
Back rest movement

Uzama Parçası
Extension Peace



Teknik Özellikler Technical Specifications

Komple Uzunluk - Full Length	1800 mm
En - Width	670 mm
Yükseklik Oturak - Seat Height	920 mm
Oturak Uzunluk - Seat Length	1140 mm
Ayak Uzatma - Foot Extension	370 mm
Sünger Dansite - Foam Density	28
Sirt Açısı - Back Angle	30°
Ayak Uzatma En - Foot Extension Width	320 mm
Ürün Ağırlığı - Product Weight	82 kg
Teker Çapı - Wheel Diameter	75 mm
Maksimum Yükseklik - Max. Height	145 mm



Arşimed Medikal Kozm. San. Tic Ltd. Şti.

Mimar Sinan Mah. Bosna Cd. No:5 / 3 Çekmeköy - İstanbul

☎ : +90 216 567 8 567 (pbx) ☎ : +90 216 567 8 367 (pbx)

🌐 : www.arsimed.com.tr 📧 : info@arsimed.com.tr



EN 60601-1 EN 60601-2



ISO 13485
certified

