

## MATTRESS SPECIFICATIONS:

- The mattress is 28kg/m<sup>3</sup> density
  - It has washable antibacterial Nylon
  - Waterproof PU cover with concealed zip on 3 sides
  - Medical grade Grey sponge
  - Dimension 190(L)x90(W)x12(H) cm
  - Mattress thickness 120 mm,
  - Latex-free, non-toxic, fire-resistant and non-carcinogenic materials
  - Antibacterial, no allergic, resistant to hospital-grade disinfection products.
- \*all mattress in our products is provided with 120 mm thickness, unless the customers required different*

Code	Description	Height (mm)	Width (mm)	Length (mm)
MATTRESS	Mattress with 28 kg/m <sup>3</sup> density	120	900	1900



• Anti-microbial

• Water resistant

• Vapour permeable cover



Arşimed Medikal Kozm. San. Tic Ltd. Şti.

Mimar Sinan Mah. Bosna Cd. No:5 / 3 Çekmeköy - İstanbul

☎ : +90 216 567 8 567 (pbx) 📠 : +90 216 567 8 367 (pbx)

🌐 : www.arsimed.com.tr 📧 : info@arsimed.com.tr

ARCHI MED<sup>®</sup>  
www.arsimed.com.tr



EN 60601-1

EN 60601-2



ISO13485  
certified